



**F**  
**l**  
**Micro**  
**W**®

---

*Quality Counts!*

**Micronucleus  
Analysis  
Kit**

**MicroFlow<sup>BASIC</sup> Archive Stage 2  
Rodent Fixed Blood**

**Instruction Manual**

---

*For Research Use Only. Not for use in diagnostic procedures.*

## Table Of Contents

1. Materials Provided.....	1
2. Additional Materials Required.....	1
3. <b>First-Time Users</b> .....	1
4. Ordering Information and Technical Services .....	1
5. Introduction.....	2
5.1. The Micronucleus Test.....	2
5.2. The MicroFlow <sup>®</sup> Method.....	2
5.3. Regulatory Acceptance.....	2
6. Ship Blood Samples to Litron .....	3
7. Results.....	3
8. References .....	4
9. License Agreement and Limited Product Warranty .....	4

## 1. Materials Provided

Kit Component	Quantity <sup>a</sup>	Storage Condition
Vacuum-insulated shipping container <sup>b</sup>	1	Ambient
Bag containing icepacks <sup>c</sup>	1	-30 °C to -10 °C
Thin clear plastic bag for shipping Study Phase Plan and Sample Submission Form	1	Ambient

- Each kit provides sufficient materials for shipping up to 60 blood samples to Litron for analysis.
- The vacuum-insulated shipping container is designed to ensure that up to 60 samples stored in LTSS will remain sufficiently cold for 48 to 72 hours.
- Upon receipt, place icepacks in a freezer set to -30 °C to -10 °C.

## 2. Additional Materials Required

- Previously fixed, frozen rodent blood samples stored at -75 °C to -85 °C in Long Term Storage Solution
- Shipping forms for overnight delivery service
- Freezer set to approximately -20 °C

## 3. First-Time Users

**We strongly recommend reading the entire instruction manual before performing these procedures.**

Please do not deviate from the procedures described in this manual. It is important that these steps are followed exactly using the supplies in this kit in order to achieve reliable results. If you have questions, please contact Litron Laboratories by calling (585) 442-0930, faxing us at (585) 442-0934, or sending an email to [info@litronlabs.com](mailto:info@litronlabs.com).

*A signed Study Phase Plan is required at Litron before sample analysis is initiated. A Sample Submission Form is also required for each shipment of samples. These forms can be found on Litron's website ([www.litronlabs.com](http://www.litronlabs.com)).*

## 4. Ordering Information and Technical Services

Litron Laboratories  
 3500 Winton Place, Suite 1B  
 Rochester, New York 14623  
 Telephone: 585-442-0930  
 Order Toll Free: 877-4-LITRON (877-454-8766)  
 Fax: 585-442-0934  
 email: [info@LitronLabs.com](mailto:info@LitronLabs.com)  
 World Wide Web: [www.LitronLabs.com](http://www.LitronLabs.com)

## 5. Introduction

This kit is used for shipping mouse or rat blood samples to Litron for flow cytometric enumeration of micronucleated erythrocyte populations.

### 5.1. The Micronucleus Test

The *in vivo* micronucleus test was established as a means of analyzing chromosomal damage. The test is based on the observation that displaced chromatin, resulting from chromosomal loss or breakage, can form a secondary nucleus (micronucleus) outside the daughter nuclei of a dividing cell. Micronuclei (MN) occur spontaneously, but an elevation in the frequency of micronuclei in a population of cells can be indicative of exposure to a genotoxic agent.

Micronuclei are particularly apparent in red blood cells (erythrocytes), which otherwise lack DNA. During erythropoiesis, a hematopoietic stem cell differentiates into an erythroblast and eventually expels its nucleus to become a reticulocyte (RET). The newly formed RET is then released from the bone marrow into the circulating bloodstream, where it develops into a mature normochromatic erythrocyte (NCE). Although the main nucleus is lost during RET formation, MN may be retained in the RET cytoplasm. Peripheral blood is ideal for micronucleus analyses because samples can be obtained from an animal easily and at multiple time points.

### 5.2. The MicroFlow<sup>®</sup> Method

Litron Laboratories has developed and patented a flow cytometric method to measure micronuclei in both the RET and NCE populations. Unlike mature NCEs, immature RETs are still rich in RNA as well as certain surface proteins (e.g., transferrin receptor, also known as CD71), and can therefore be differentially stained based on these features. An increase in the frequency of micronucleated reticulocytes (MN-RETs) can indicate acute genotoxicity associated with a recent cell division. In mice, an increase in the frequency of micronuclei in the NCE population (MN-NCE) can indicate accumulated DNA damage associated with a sub-chronic or chronic treatment regimen. Elevated MN-NCE frequencies in rat blood need to be interpreted with caution, since splenic filtration function is the dominant factor that influences these values.

The MicroFlow method offers significant advantages compared to traditional microscopic scoring, such as:

- Greater number of cells can be examined for MN
- Faster data acquisition
- Increased statistical power of the assay
- Objective analysis of samples

The MicroFlow method also offers advantages over other automated methods, including:

- Availability for many species of toxicological interest
- Anti-platelet antibody to ensure reliable data
- Calibration Standards to ensure intra- and inter-laboratory reproducibility of data
- Ability to store samples for extended periods of time before analysis

Crucial components of this method are the Calibration Standards, which aid flow cytometer configuration for the micronucleus scoring application. Fixed blood from animals infected with *Plasmodium berghei* are used to configure the flow cytometer before analysis. Whereas MN are relatively rare and exhibit a heterogeneous DNA content, parasitized cells are prevalent and have a homogenous DNA content. These characteristics make them ideal for calibrating the flow cytometer for the micronucleus scoring application. After optimizing the flow cytometer with the Calibration Standards, micronucleus analyses can be performed reliably and with minimal intra- and inter-experimental variation.

### 5.3. Regulatory Acceptance

The US FDA accepts preclinical MicroFlow data, and this method adheres to the necessary guidelines as stated by the International Workshop on Genotoxicity Test Procedures (IWGTP). Additionally, the most current Organization for Economic Co-Operation and Development (OECD) guidelines regarding micronucleus testing, Guideline 474, indicates that flow cytometry, using appropriate calibration standards, can provide better inter- and intra-laboratory reproducibility and sensitivity than manual microscopic scoring. It also states that “Commonly used laboratory strains of healthy young adult animals should be employed. Mice, rats, or another appropriate mammalian species may be used. When peripheral blood is used, it must be established that splenic removal of micronucleated cells from the circulation does not compromise the detection of induced micronuclei in the species selected. This has been clearly demonstrated for mouse and rat peripheral blood.”

## 6. Ship Blood Samples to Litron

1. Before preparing shipment, verify that all icepacks have been frozen in a freezer set to approximately -20 °C. Please use the shipping container that was provided by Litron. It is imperative that all icepacks are frozen and that the Litron-provided vacuum-insulated shipping container is used in order to maintain correct temperatures during transit.
2. **Complete Study Paperwork and Prepare Shipment.** Complete and sign the appropriate documents and place them inside the thin clear plastic bag. Sample analysis cannot begin until the Sample Submission Form and Study Phase Plan are received.

Position the following items close to the freezer where samples are stored:

- Vacuum-insulated shipping container (supplied with kit)
- Shipping bag (supplied with kit)
- Frozen Ice packs (supplied with kit)

3. **Place the bagged cryovial storage boxes upright at the bottom of the shipping box** (up to four boxes should fit in one layer on the bottom of the box). Do not ship more than eight cryovial storage boxes in any one box. Fill with frozen ice packs. Place the cardboard piece, silver panel, and foam piece on top of the icepacks (in this order). If you need additional shipping boxes or icepacks, contact Litron.
4. **Seal Box.** Place the thin clear plastic bag containing the applicable forms on top of the foam piece. Close the cardboard flaps of the outer box and use shipping tape to secure the box top. *It is very important to get a tight seal in order to maintain temperatures during transit.*
5. **Ship To Address.** *DO NOT Identify the samples as Category A or B biological substances.* Please verify that the shipper you use guarantees overnight delivery, and ship to the following address:  
Litron Laboratories  
3500 Winton Place, Suite 1B  
Rochester, New York 14623  
585-442-0930

Unexpected shipping delays may occur at any time. Therefore, it is best to ship samples on Monday or Tuesday and avoid shipping during holidays.
--

## 7. Results

Preliminary results will be emailed and a hard copy of the final results will be provided, if requested.

## 8. References

- Asanami S, Shimono K, Sawamoto O, Kurisu K and Uejima M (1995) *Mutation Research* 347, 73-78.
- Asano N, Torous DK, Tometsko CR, Dertinger SD, Morita T, and Hayashi M (2006) *Mutagenesis* 21(1), 15-20
- Wakata A, Miyamae Y, Sato S, Suzuki T, Morita T, Asano N, Awogi T, Kondo K and Hayashi M (1998) *Environmental Molecular Mutagenesis* 32, 84-100.
- Dertinger SD, Torous DK and Tometsko KR (1996) *Mutation Research* 371, 283–292.
- Dertinger S, Torous D and Tometsko K (January 1999 [filed September 1996]) U.S. Patent No. 5,858,667.
- Dertinger SD, Torous DK and Tometsko KR (August 2000 [filed January 1999]) U.S. Patent No. 6,100,038.
- Dertinger, S (Patent Pending [filed 6/28/04]).
- Dertinger SD, Torous DK, Hall NE, Tometsko CR and Gasiewicz TA (2000) *Mutation Research* 464, 195–200.
- Dertinger SD, Camphausen K, Macgregor JT, Bishop ME, Torous DK, Avlasevich S, Cairns S, Tometsko CR, Menard C, Muanza T, Chen Y, Miller RK, Cederbrant K, Sandelin K, Pontén I and Bolcsfoldi G (2004) *Environmental Molecular Mutagenesis* 44, 427-435.
- Dertinger SD, Bishop ME, McNamee JP, Hayashi M, Suzuki T, Asano N, Nakajima M, Saito J, Moore M, Torous DK and MacGregor JT (2006) *Toxicological Sciences* 94, 83-91.
- Dertinger SD, Tsai Y, Nowak I, Hyrien O, Sun H, Bemis JC, Torous DK, Keng P, Palis J and Chen Y (2007) *Mutation Research* 634, 119-125.
- Hayashi M, Sofuni T and Ishidate M (1983) *Mutation Research* 121, 241-247.
- Hayashi M, MacGregor JT, Gatehouse DG, Blakey DH, Dertinger SD, Abramsson-Zetterberg L, Krishna G, Morita T, Russo A, Asano N, Suzuki H, Ohyama W and Gibson D (2007) *Mutation Research* 627, 10-30.
- Heddle JA (1973) *Mutation Research* 18, 187–190.
- Heddle JA, Hite M, Kirkhart B, Mavournin K, MacGregor JT, Newell GW and Salamone MF (1983) *Mutation Research* 123, 61-118.
- Holden H, Majeska J and Studwell D (1997) *Mutation Research* 391, 87-89.
- Kissling GE, Dertinger SD, Hayashi M and MacGregor JT (2007) *Mutation Research* 634, 235-240.
- MacGregor JT, Wehr CM and Gould DH, (1980) *Environmental Mutagenesis* 2, 509–514.
- MacGregor JT, Bishop ME, McNamee JP, Hayashi M, Asano N, Wakata A, Nakajima M, Saito J, Aidoo A, Moore M and Dertinger SD (2006) *Toxicological Sciences* 94, 94-107.
- Salamone M, Heddle J, Stuart E and Katz M (1980) *Mutation Research* 74, 347-356.
- Salamone MF and Heddle JA (1983) In: F.J. de Serres, ed. *Chemical Mutagens: Principles and Methods for their Detection*, Vol 8. New York: Plenum, 1983; 111-149.
- Schmid W (1975) *Mutation Research* 31, 9-15.
- Serke S and Huhn D (1992) *British Journal of Haematology* 81, 432-439.
- Tometsko AM (July 1993 [filed March 1990]) U.S. Patent No. 5,229,265.
- Tometsko AM, Torous DK and Dertinger SD (1993) *Mutation Research* 292, 129–135.
- Tometsko AM, Dertinger SD and Torous DK (1993) *Mutation Research* 292, 137-143.
- Tometsko AM, Torous DK and Dertinger SD (1993) *Mutation Research* 292, 145-153.
- Tometsko AM, Dertinger SD and Torous DK (1995) *Mutation Research* 334, 9-18.
- Torous DK, Hall NE, Dertinger SD, Diehl MS, Illi-Love AH, Cederbrant K, Sandelin K, Bolcsfoldi B, Ferguson LR, Pearson A, Majeska JB, Tarca JP, Hewish DR, Doughty L, Fenech M, Weaver JL, Broud DD, Gatehouse DG, Hynes GM, Kwanyuen P, McLean J, McNamee JP, Parenteau M, Van Hoof V, Vanparys P, Lenarczyk M, Siennicka J, Litwinska B, Slowikowska MG, Harbach PR, Johnson CW, Zhao S, Aaron CS, Lynch AM, Marshall IC, Rodgers B and Tometsko CR (2001) *Environmental and Molecular Mutagenesis* 38, 59–68.
- Torous DK, Hall NE, Murante FG, Gleason SE, Tometsko CR and Dertinger SD (2003) *Toxicological Sciences* 74, 309-314.
- Torous DK, Hall NE, Illi-Love AH, Diehl MD, Cederbrant K, Sandelin K, Pontén I, Bolcsfoldi G, Ferguson LR, Pearson A, Majeska JB, Tarca JP, Hynes GM, Lynch AM, McNamee JP, Bellier PV, Parenteau M, Blakey D, Bayley J, van der Leede BM, Vanparys P, Harbach PR, Zhao S, Filipunas AL, Johnson CW, Tometsko CR and Dertinger SD (2005) *Environmental Molecular Mutagenesis* 45, 44-55.
- Torous D, Asano N, Tometsko C, Sugunan S, Dertinger S, Morita T and Hayashi M (2006) *Mutagenesis* 21, 11-13
- Torous DK, Dertinger SD, Hall NE and Tometsko CR (2000) *Mutation Research* 465, 91-99.
- Witt KL, Livanos E, Kissling GE, Torous DK, Caspary W, Tice RR, and Recio L (2007) *Mutation Research* 649, 101-113.

## 9. License Agreement and Limited Product Warranty

By utilizing this kit, your company is agreeing to be bound by the terms of this License. This License allows the use of the MicroFlow<sup>®</sup> Kit for the analysis of 60 samples, either in-house (MicroFlow<sup>PLUS</sup> Kit), or at Litron's facility (MicroFlow<sup>BASIC</sup> Kit).

MicroFlow<sup>®</sup>. All rights reserved. *MicroFlow*<sup>®</sup> is a registered trademark of Litron Laboratories. Patent Nos. 7,425,421, 7,867,447, 8,076,095, 8,586,321, 8,889,369, 9,285,365, 2,529,802, and 1,649,038. Copyright 2003-2017, Litron.

By accepting these products, you acknowledge that they will be used in accordance with their intended labeling (For *in vitro* research use only. Not for human or animal diagnostic or therapeutic use.). Uses other than the labeled intended use may be a violation of local laws.

This warranty limits our liability to replacement of this product. Litron shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.